

Vacuum drying chambers



Product parameters

Brand	INNOVA	INNOVA	INNOVA	INNOVA
Model	INO-25/ZK	INO-35/ZK	INO-45/ZK	INO-55/ZK
Inner Chamber Size	275x300x300	345x415x370	450x450x450	550x650x600
Liner Material	Stainless steels	Stainless steels	Stainless steels	Stainless steels
Voltage	220	220	220	220
Power	0.9	1.4	1.6	2
Temp. Range	RT+10~200	RT+10~200	RT+10~200	RT+50~200
Accuracy	±1	±1	±1	±1
Stable Volatility	±2	±2	±2	±2
Volatility	±2	±2	±2	±2
Vacuum	133Pa	133Pa	133Pa	133Pa
Ex. grade	ExdIIBT4Gb	ExdIIBT4Gb	ExdIIBT4Gb	ExdIIBT4Gb