

ULT FREEZER & MEDICAL FRIDGE



Qingdao Innova Bio-Meditech Co., Ltd.

Tel.: +86 532 8789 0634/Email: info@innobiomed.com

Product Feature

Temperature Control	▶	Microcomputer digital control. Temperature realization display. Display accuracy is 0.1 °C
	▶	The temperature can be adjusted from -40 to -86 degrees.
	▶	With 7-inch LCD touch screen, easy to operate, easy to observe, clear and intuitive;
	▶	Alarm Temperature can be set.
Security System	▶	With a perfect alarm system. Showing with sound and light .
	▶	Display input voltage, with USB output, support data archive / 1 month, data view, with 485 interface, can be connected to a printer
	▶	High and low temperature alarm, sensor failure alarm, power failure alarm, door opening alarm, low battery, abnormal ring temperature alarm and other functions can be realized
Refrigeration System	▶	The SECOP compressor make the performance stable.
	▶	The condensing fan brand in EBM, which made the noise in low.
	▶	Ultra-thick foam layer, double sealing of the door body, effective protection of the temperature inside the freezer.
	▶	100% CFC free.
Ergonomic design	▶	Easy Access Built-in 304 stainless steel shelves can be adjusted.
	▶	High Security The door is equipped with a security lock.
	▶	High Applicability 187V-242V wide voltage design, suitable for unstable voltage environment.
Certifications	▶	CE, FDA, ISO9001, ISO13485

Product Parameters

Item No.	MF-86V508T	Rated Input(w)	880
Capacity/L	508	Power Consumption	12.9 KWh/24H
Temperature	-40~-86℃	Thermostat	Digital & LCD Touch Screen
Power Supply/V/Hz	220/50	Display	
Refrigerant	mixed with CFC free	Fan of Condenser	EBM
Noise /dB(A)	60	Cooling Way	Direct Cooling
Compressor Brand	SECOP	Inner Size/mm	570*710*1225
Compressor Quantity	2	Outer Size/mm	870*1030*1975
Shelves	3	Packing Size/mm	950*1130*2125
Doors	4	Foaming thickness/mm	150
Lock	YES	Inner Material	304 stainless steel
Wheels	YES	Weight/kgs	250
Alarm Way	Sound+Light	Defrosting Way	Manual
Alarm Type	High and low temperature alarm, sensor failure alarm, power failure alarm, door opening alarm, low battery, abnormal ring temperature alarm		
Packing	Plywood	Warranty	1 year for machina 3 years for compressor
USB	YES	LCD Screen	YES

Product Pictures

